

Notification

Fixed Vulnerabilities When Image Reading in NIS-Elements Viewer.

Published on August 29, 2024

Healthcare Business Unit

NIKON CORPORATION

1. Overview

Vulnerabilities has been Fixed in NIS-Elements Viewer when reading image files.

2. How to check if your software affected:

The affected products and version are as follows:

| Product name | Version |
|---------------------|-------------------------|
| NIS-Elements Viewer | Ver.5.21.00 and earlier |

To check your software's version, refer to the steps below:

Go to Menu and select [Help] > [About], The Version information dialog box will open.

3. Fixed Vulnerabilities

The following vulnerabilities has been fixed.

| Vulnerability Category | Vulnerability Impact | CVE-ID | Affected Image Files |
|----------------------------|------------------------|--------------------------------|----------------------|
| out-of-bounds read | execute malicious code | CVE-2022-40656 | ND2 file |
| | | CVE-2022-40662 | TIFF file |
| | | CVE-2022-40663 | TIFF file |
| heap-based buffer overflow | execute malicious code | CVE-2022-40661 | BMP file |
| | | CVE-2022-40660 | PSD file |
| | | CVE-2022-40655 | ND2 file |
| out-of-bounds write | execute malicious code | CVE-2022-40657 | PSD file |
| | | CVE-2022-40658 | TIFF file |
| | | CVE-2022-40659 | TIFF file |

4. Threats posed by the vulnerabilities

Opening an attacker-crafted image files may allow arbitrary code to be executed.

5. Countermeasure

The countermeasure is as follows:

| Product name | Version |
|---------------------|---|
| NIS-Elements Viewer | Update the software to Ver.5.22 or later. |

Please download the software from the following URL:

[Europe and Asia Region]

https://www.microscope.healthcare.nikon.com/en_AOM/software-firmware

[Americas Region]

<https://www.microscope.healthcare.nikon.com/products/software/nis-elements/viewer>

6. How to reduce/avoid the threat

If you are unable to upgrade your product immediately, we recommend the following mitigation measures:

- Do not open untrusted files.

7. Contact information

Please contact your sales representative.

The information is current as of the date of publication. It is subject to change without notice.