Notification

Addressing PHP Vulnerabilities in NIS-Elements StatViewer

Published: October 28, 2026

Healthcare Business Unit
NIKON CORPORATION

1. Overview

Security vulnerabilities have been identified in the PHP component used by NIS-Elements StatViewer, which is bundled with NIS-Elements. No adverse effects on our products have been confirmed at this time. However, customers using the affected products are requested to take the actions described below.

2. Affected Products and Versions

The following versions of NIS-Elements include NIS-Elements StatViewer.

Product Name	Version
NIS-Elements Advanced Research	Ver.4.00–Ver.6.10.03
NIS-Elements Basic Research	Ver.4.00–Ver.6.10.03
NIS-Elements Documentation	Ver.4.00–Ver.6.10.03
NIS-Elements Confocal	Ver.4.00–Ver.6.10.03
NIS-Elements Enhanced Resolution	Ver.4.00–Ver.6.10.03

Note: NIS-Elements StatViewer is not bundled with NIS-Elements Ver.6.20 or later.

3. How to Check

You can check whether NIS-Elements StatViewer is installed as follows.

1. In Windows, open Settings → Apps → Installed apps.
2. Look for “NIS-Elements StatViewer *.*” in the list or use the search box.

4. Remediation

If NIS-Elements StatViewer is installed on your system, please uninstall (remove) it.

Uninstalling NIS-Elements StatViewer does not affect the functionality or operation of NIS-Elements itself.

4.1 Uninstalling Using Standard Windows Features

1. In Windows, open Settings → Apps → Installed apps. Depending on your Windows version, open Apps & features, or Programs and Features in Control Panel.
2. Look for “NIS-Elements StatViewer *.*”.
3. Select Uninstall and follow the on-screen instructions.

4.2 Manual Removal

If you cannot uninstall the software using the method in Section 4.1, you can remove it manually by following the steps below. Administrator privileges are required for manual removal.

1. Run Command Prompt as an administrator.
2. Run the following commands to stop and delete the StatViewer service.

```
net.exe stop nis-e_stats
sc.exe delete nis-e_stats
```

3. Close Command Prompt.
4. Delete the StatViewer installation folder. Typical installation locations are listed below.

- C:\Program Files\NIS-Elements StatViewer
- C:\Program Files\NIS-Elements D\StatViewer
- C:\Program Files\NIS-Elements\StatViewer

Note: If the software is installed in a location other than those listed above, check the applicable installation folder.

5. Affected Vulnerabilities

The vulnerabilities associated with the PHP component used by NIS-Elements StatViewer are listed below.

Vulnerability	CVE ID	Affected Software
Vulnerabilities related to PHP2Shell	CVE-2026-109618 CVE-2026-109619	NIS-Elements StatViewer

Under certain conditions, these vulnerabilities could allow arbitrary code to be executed remotely without authentication. If exploited, they could lead to control of the server or information disclosure.

6. Contact Information

Please contact your sales representative.

The information posted here is current as of the date of announcement. Please note that the information may differ from the latest information, such as products that are no longer on sale or organizational changes.